

FootCare





URIEL®

SUPPORTS



We all need
a little support

Wir alle benötigen ein wenig Stütze
Todos necesitamos un poco de apoyo
Nous necessitions tous un peu de support




FootCare

URIEL
Sport

HTS

PRESSURE
GARMENT

MEASUREMENTS



INTRODUCTION - COMPANY PROFILE



"URIEL" stands for high-quality health supports. "MEDITEX Limited" was founded in 1975 by Mr. Uriel Rosen O.M. in Israel. He received his training in Germany. Today, MEDITEX Ltd. sells its products in many countries around the world and is proud of its wide range of products.

This catalog presents the complete "URIEL" line of Health Supports and related products and is suitable for the retail market, as well as for the orthopedic institutes and medical offices. All products are manufactured under the control of Mr. Rony Rosen (son of U. Rosen) and with the experienced guidance of medical experts. MEDITEX offers outstanding and consistent quality, backed by effective and on going product research to develop and renew our products constantly. Superior quality materials, expert manufacturing techniques and innovative product designs have earned "URIEL" an international reputation. Advances in knitting technology and elastic yarn formulations enable the "URIEL" range to continuously bring orthopaedically correct, therapeutically designed devices to the market.

This catalog will guide you to choose the most suitable device for the injured or handicapped patients. It will familiarize you completely with the "URIEL" Health Supports and will further aid you in your professional work. Our recommendations may assist you to select the right products for each case.



Uriel Rosen O.M.
General manager



Rony Rosen
Vice General manager
Head of Pressure Garment Dept.

Order CODE SYSTEM

PACKING	MATERIAL	ITEM No.	EXPLA.	COLOUR	SIZE
☐ = BOX ☒ = POLY BAG	☐ = MEDICAL (knitted) ☒ = NEOPRENE ST = SPORT (Neoprene) H = H.T.S. N = SUPORTO F = NEO FIX	T 4 3	M = MEN W = WOMEN L = LEFT SIDE R = RIGHT SIDE U = UNIVERSAL A = MACH MATERIAL E = LONG / EXTENDED H = MONO CENTRIC D = DUAL CENTRIC G = DUAL CEN. ADJUST. P = PELOTA / INSERT / HORSESHOE	B = BLACK BL = BLUE C = SKIN G = GRAY R = RED W = WHITE Y = YELLOW	B = BLACK BL = BLUE C = SKIN G = GRAY R = RED W = WHITE Y = YELLOW



CERVICAL FOAM
CERVICAL FIRM
SLEEP PILLOW

COLLAR

PAGE 3



ELBOW SUPPORT
TENNIS BANDAGE
ARM SLING POUCH STRAP
ARM+SHOULDER IMMOBILIZER

ELBOW

PAGE 4-7



WRIST SPLINT
WRIST SUPPORT
THUMB SPLINT

WRIST


PAGE 8-10



CLAVICULA BANDAGE
SHOULDER SUPPORT

SHOULDER

PAGE 11




MATERNITY BELT
INDUSTRIAL BELT
ABDOMINAL BELT

MAGNETO BELT
SURGICAL BELT
LIGHT MEN'S BELT

LUMBO BELT
WOMEN'S BELT
THERMO BELT

BELT

PAGE 12-17



ABDUCTION BANDAGE
PAVLIK BANDAGE

THIGH CAP
SUSPENSORY
HERNIA BELT

HIP BANDAGE

PAGE 18-19



HIP PROTECTOR

KNEE BRACE
KNEE SUPPORT
KNEE IMMOBILIZER
HINGED KNEE BRACE

KNEE

PAGE 20-25



ANKLE BRACE
ANKLE SUPPORT
ANKLE AIR BRACE
WALKERS FOOT WEAR

ANKLE

PAGE 26-29



INSOLES
ARCH SUPPORT
HALLUX VALGUS SPLINT
SILICONE HEEL CUSHIONS

PEDICURE PRODUCTST
DIBETEIC SILVER SOCKS
SILICONE FOOTCARE PRODUCTS

FOOT CARE

PAGE 30-37



PRESSURE GARMENTLIPOSUCTION

PAGE 38-39

BANDAGE

PAGE 40

NEOPRENE BANDAGE ONE SIZE FITS ALL

PAGE 41

**PRESSURE
GARMENT**

HIS

**ORIEL
Sport**

MEASUREMENTS

PAGE 42-43

MEASUREMENTS



MEDICAL SUPPORT



HEALTH SUPPORTS

WE ALL NEED
a little
SUPPORT



SUPPORTO

The range of sport medicine compression supports and general supporting appliances are based on knitted, sweat-absorbing natural cotton, non-irritating to the skin.

See page 30



LIPOSUCTION
SPECIAL CORSETS, AFTER PLASTIC SURGERY



See page 38



SPECIAL DEPARTMENT FOR THE PRODUCTION OF
PRESSURE GARMENT



See page 48




**THERMO LINE
NEOPRENE SUPPORT**

HEAT RETAINING MATERIAL

The Uriel support braces are manufactured from 5 mm NEOPRENE rubber with a nylon laminate and 100% cotton terry cloth lining, based upon a heat retaining treatment. NEOPRENE insulates the

body's natural warmth, aids blood circulation, helps to reduce edema and allows the soft tissues to become less susceptible to injuries, although severe injuries may need stiff braces.

**SYSTEM OF
NEOPRENE
HEAT
RETAINER**



See page 41

**Neoprene
bandage**



Sport

ONE SIZE FITS ALL

UNE TAILLE POUR TOUS

EINE GRÖSSE - PASST ALLEN

TALLE PARA TODAS LAS MEDIDAS



See page 30



100% SILICONE



HI-TEX-SUPPORT IS MADE OF AN INNOVATIVE, ESPECIALLY THIN **LYCRA** FABRIC. OPTIMAL FITTING THE INDIVIDUAL ANATOMICAL SHAPE OF THE BODY LIKE A SECOND SKIN ACCORDING TO THE MEMORY EFFECT

SYSTEM OF



MATERIAL COMPOSITION:

42% **LYCRA**
16% **POLYURITHAN**
42% **POLYMICRO FLEECE**



LYCRA FABRIC
POLYURITHAN
POLY FLEECE

BREATHABLE

See page 32



60/62 - FOAM CERVICAL COLLAR

1. Light neck injury.
2. Suitable for use during sleep.



ITEM 62
ANATOMICAL CONTOUR COLLAR



60

62



60/62/65

Available in 100% cotton stretch cover
Available in 100% cotton towel cover

colors: white, skin,
colors: white, skin, grey and yellow

COLOR

W

C

G

Y

65 - FIRM CERVICAL COLLAR

1. Suitable for mild to chronic cases of neck injury.
2. Firm when applied, yet pliable enough to be comfortable.

The patient should be able to see his/her feet when the collar is applied.



65



63 - PLASTIC COLLAR

Adjustable height extension

Suitable for mild to chronic cases of neck injury.

The patient must always have the ability to see his/her feet when the collar is applied.



63



66 - PHILADELPHIA COLLAR

This emergency collar is indispensable in first-aid-treatment, in order to prevent irreversible injuries, since it provides an excellent cervical spine immobilization and also offers adequate access to the mouth for resuscitation and airway maintenance procedures such as suctioning. Can be easily applied to patients in either sitting or supine position.

66






Measurements on page 42

67 - DREAM PILLOW

1. Very helpful for painful necks.
2. Its wavy form helps to ease head positions.
3. The cavities at the back prevent pressure on the ears while sleeping on the side, preventing the need to put the hand under the ear, which may cause blood circulation disturbances.

67



MEASUREMENTS



4



URIEL[®] ARM SLING



FootCare

URIEL[®]
Sport

HIS

PRESSURE
GARMENT

MEASUREMENTS

85



85 - POUCH ARM SLING

1. Holds the arm securely.
2. Light material designed for warm climates.
3. Skin friendly.
4. Adjustable strap and buckle.

85+



85+ - FIXATION ARM SLING Pouch

As item 85 but with additional
strap for arm immobilization.

80



80 - ARM SLING STRAP

1. Holds the arm securely.
2. Fits right or left arm.
3. Can be worn with or without cast.
4. Adjustable strap and buckle.

82



82 - TUBE ARM SLING

(Adjustable)

1. Holds the arm securely.
2. Extra comfortable sling.
3. Polyfoam and cotton cover.
4. Minimal pressure and tension on the back of the neck.
5. Velcro fasteners for easy application.

SHOULDER SUPPORT

T78 - SHOULDER SUPPORT

T78

For relief of numerous injuries,
Rotator Syndrom, Bicipital Tenditis, Acromio
Clavicular, Osteoarthritis.
Fits left and right arm, made of Neoprene laminate.

Made of :
NEOPRENE LAMINATE



86 - ARM SHOULDER BANDAGE

86

1. Limits arm and shoulder movement.
2. Adjustable with velcro fastener.

Fits left and right arm



87 - ARM & SHOULDER IMMOBILIZER

87

- Used post-operative.
1. Immobilizes the shoulder and upper arm region.
 2. Adjustable with Velcro fastener

Fits left and right arm



FootCare

URIEL
Sport

HIS

PRESSURE
GARMENT



**PRESSURE
GARMENT**

6



TENNIS BANDAGE



T90



T90 - TENNIS BANDAGE SUPPORT

1. Relief from tennis elbow condition.
2. Absorbent terry cloth.
3. Pressure on common extensor.



with plastic insert,
pressure on extens
or muscles.

made of : neoprene laminate



91

91 - TENNIS BANDAGE

1. Relief from tennis elbow condition.
2. Absorbent terry cloth.
3. Pressure on common extensor.

911



with silicone cushion

911 - TENNIS BANDAGE

As item 91

1. Relief from tennis elbow condition.
2. Made of absorbent terry fabric.
3. The silicone cushion exerts compression on common extensor.
4. Inserted elastic strap allows pressure adjustment.

One size fits all.



95+

95+ - ELBOW BANDAGE

1. Relief from tennis elbow condition.
2. Pressure on common extensor
3. as item 95 with tennis elbow strap.

Fits left and right arm.



T95+

T95+ - ELBOW SUPPORT

1. Relief from tennis elbow condition.
2. Pressure on common extensor.
3. As item T95 with tennis elbow strap.

Made of :
NEOPRENE LAMINATE

Fits left and right arm.

Measurements on page 42

MEASUREMENTS



ELBOW BANDAGE



7



**PRESSURE
GARMENT**

922 - ELBOWSIL BANDAGE

1. Relief from tennis elbow condition.
2. Pressure on common extensor.
With 2 silicone pads.
3. For relief of painful elbow joints caused by overextension, strain and fatigue.
4. Compression to arthritic joints.
5. Anatomically contoured, with silicone cushions along the radial and ulnar side.

Fits left and right arm.



95 - ELBOW BANDAGE

1. For relief of painful elbow joints caused by overextension, strain and fatigue.
2. Heat retention.
3. Compression to arthritic joints.

95



T95 - ELBOW SUPPORT

1. For relief of painful elbow joints due to overextension, strain and fatigue.
2. When applied, the support provides warmth and compression to arthritic sufferers.

MADE OF :
NEOPRENE LAMINATE

T95



T98 - ELBOW HYPER-EX SUPPORT

1. Prevents Hyper Extension of the elbow joint.
2. Use after olecranon fracture.
3. For use during sports activity.

Made of :
NEOPRENE LAMINATE

T98



Measurements on page 42

MEASUREMENTS



8



WRIST BANDAGE



**PRESSURE
GARMENT**



20

20 - WRIST BANDAGE

1. Prevents wrist aches.
2. Post injury wrist support.
3. Light material designed for warm climates.
4. Hook and loop fasteners for easy application.
5. Adjustable sizes due to fasteners.

25 - WRIST SUPPORT



25

T25 - WRIST SUPPORT



T25

1. Adjustable thumb strap, helps to support the base of the thumb.
2. The support can be worn in a number of positions due to the stabilizer - strap.
3. Easier to apply than conventional wrist support with 3 tabs for closure of the support.
4. Prevents wrist aches.

Made of :
NEOPRENE LAMINATE



21

21 - WRIST BANDAGE

1. Supports the wrist, gives more tension and adjustment due to an additional strap.
2. Prevents wrist aches.
3. Use against swellings.



NEW



222R

222 - WRIST SILICONE BANDAGE

222L

222R RIGHT
222L LEFT

with silicone cushion.

1. Supports the wrist, gives more tension and adjustment due to an additional strap.
2. Prevents wrist aches.
3. Use against swellings.
4. With silicone cushion.

WRIST SPLINT

26 - WRIST SPLINT

SUITABLE FOR BOTH SIDES

Includes hardened aluminium splint, designed as cock-up splint.



This splint is adaptable for right and left hand use. To alter sides, remove the insert, rotate it 180 degrees and re-insert.

1. Recommended in cases of fractures, after removal of a cast or injuries to ligaments, muscles and tendons.
2. It is designed to provide a strong support and to minimize movement of the wrist.
3. The splint holds the hand at an angle of approximately 35 degrees to the wrist. The angle can be altered by slightly bending the aluminium splint.

T26

Made of :
NEOPRENE LAMINATE



26



2722 R/L - WRIST SPLINT WITH SILICONE CUSHION

2722 L

2722 R

NEW

As item 26, with stronger fixation due to an additional fastener strap, a silicone pad and an additional support-rod at the back of the hand

Made of:
bi-elastic knitted material



with silicone cushion.



27 - WRIST SPLINT

27L - Left 27R - Right

As item 26. Recommended for the physically handicapped who have difficulties closing the fastener.

Includes hardened aluminium splint, designed as cock up splint.

27 L

27 R



T 27 L

T 27 R



NEOPRENE LAMINATE

27E - LONG WRIST SPLINT

27EL - Left 27ER - Right

As item 27 but longer and with additional fixation strap.

Includes hardened aluminium splint, designed as cock up splint.

27E L

27E R



MEASUREMENTS



10 

URIEL

THUMB SPLINT



278

278 - THUMB AND WRIST SPLINT

This splint is a combination of item 27 and item 28.

It restricts the movement of the thumb and supports the wrist due to the inserted aluminum rod.

Useful for the treatment of sprained or strained wrists and thumbs and of inflammation in this area.

With adjustable thumb strap.

Fits either left or right hand.



28

28 - T28 THUMB SPLINT

1. This splint allows minimal movement of the thumb and supports the wrist, allowing finger movement.

2. Useful for the treatment of sprained or strained thumbs & wrists.

3. Indications: Trigger Thumbs-M.P. joint Carpo First metacarpal joint, Skier's thumb, Game-keepers thumb, Bowlers thumb, Tenosynovitis, non-specific thumb or wrist pain.

4. This splint can be used as a resting splint after injury or surgery. Local compression and support is achieved around the wrist and thumb area.

5. Fits either left or right hand.

6. Useful for the treatment of inflammation of the wrist and thumb area.

Made of : NEOPRENE

T28



288

T288 - T288 THUMB SPLINT

With adjustable Thumb strap

NEW

Made of : NEOPRENE

T288



T28A

T28A - THUMB SPLINT

1. This splint allows minimal movement of the thumb.

2. Useful for the treatment of sprained or strained thumbs.

3. Indications: Trigger Thumbs-M.P. joint Carpo First metacarpal joint, Skier's thumb, Game-keepers thumb, Bowlers thumb, Tenosynovitis, non-specific thumb pain.

4. This splint can be used as a resting splint after Injury or surgery. Local compression and support is achieved around thumb area.

5. Fits either left or right hand.

Made of :
NEOPRENE LAMINATE

FootCare

URIEL
Sport
HIS

PRESSURE
GARMENT

MEASUREMENTS

Measurements on page 42

74 - CLAVICLE SUPPORT

This figure "8" strap, ensures a convenient conformity in the axilla, alleviates pressure on the brachial artery or nerve.

Available in 2 sizes S/M - L/XL


74



76 - SHOULDER SUPPORT

1. Improves posture while allowing freedom of body movement.
2. Suitable to be worn directly on the skin due to special non-skin irritating material.
3. Can be used as psychological and physiological support for proper posturing during adolescence.
4. Light material designed for warmer climates.
5. Material used is very comfortable on the skin.
6. Hook and loop fasteners for easy application.
7. Adjustable sizes due to fasteners.
8. Used for the treatment of round shoulders in all stages of life.

75



76



776 - LUMBO DORSAL BELT

Specially for older women with problems caused by osteoporosis.
Support for upper and lower back.

776



**PRESSURE
GARMENT**



URIEL
Sport



PRESSURE
GARMENT

12

URIEL® ABDOMINAL BELT



11

11 - ABDOMINAL AIR BELT

(Width: 20 cm / 8")

1. Supports the abdomen - post surgical application.
2. Useful in cases of sagging abdomen.
3. Improves body shape.
4. Useful as a post pregnancy brace.
5. Hook and loop fasteners for easy application & adjustable sizes.

Ventilated, knitted material.



T11C

T11C - ABDOMINAL BELT

Width: 20 cm/10"

1. Supports the abdomen - post surgical application.
2. Useful in cases of sagging abdomen.
3. Improves body shape.
4. Useful as a post pregnancy brace.
5. Hook and loop fasteners for easy application & adjustable sizes.

Made of : NEOPRENE LAMINATE



11EU

11EU - ABDOMINAL BELT

(Width: 25cm / 10")

1. Supports the abdomen - post surgical application.
2. Useful in cases of sagging abdomen.
3. Improves body shape.
4. Useful as a post pregnancy brace.
5. Hook and loop fasteners for easy application & adjustable sizes.

Made of: POLYESTER.



11PS

11PS - POST SURGICAL BELT (Width: 30 cm / 12")

1. Specifically designed for post surgical abdominal application.
2. **The material is 100% cotton,** allowing use immediately after accidents in cases of abdominal injuries.
3. Extremely comfortable material, to wear on the skin.
4. Keeps stitches and pads in place.
5. Suitable for open wound dressing.

04U - LUMBO BELT (light)

(Made of **light** material)

Back part 15 cm high.

Relieves pain in the back caused by over-stress to the coccyx.

Useful for treatment of pain in the lower back.

04U



06WF - WOMEN LUMBAR BELT - FORTE

Specially designed for women

Made of **strong** material

Back part 25 cm high.

designed for wider hips and smaller waist.

06W - WOMEN LUMBAR BELT - LIGHT

ES ITEM 06WF

Made of **light** material

Back part 25 cm high.

06W
light



06WF



0622 - LUMBOSIL SACRAL BELT

1. Useful for treatment of pain in the lower back.

2. With silicone pelota.

Made of:
bi-elastic knitted material

NEW

with silicone cushion

0622



07AE - DORSO SACRAL (DLS)B.

Back part 30 cm / 12" high.

1. Useful for treatment of pain in the lower back.

2. Made of bi-elastic, knitted cotton light porous material.

3. 4 support rods situated along the spine.

07AE



19

STRONG

MATERNITY BELT

1. Lifts the abdomen to prevent lower back pain.
2. Includes a strap that aids in the lifting process of the abdomen.
3. Recommended for use from the 6th month of pregnancy.
4. The wearing of this brace does not disturb the growth process of the fetus.



18

LIGHT



PRESSURE GARMENT



MEASUREMENTS

14 



LIGHT BACK BELT



10

10 - LUMBO BELT

(Width: 20cm / 10")

1. Useful for treatment of pain in the lower back.
2. Made of bi-elastic, knitted cotton mat.
3. Equipped with 4 support rods situated along the spine.
4. Belt alleviates tension on the spine by controlling abdominal sag.



06

06 - MICDO BELT


(Width: 20cm / 10")

(Made of light material)

Back part 25 cm high.

Support for the lower back.

1. Useful for treatment of pain in the lower back.
2. Equipped with 4 support rods situated at back of spine.



12

12 - MAGNETO BELT

(Width: 20cm / 10")

Back part 25 cm high.

1. Treatment of mechanical low back pain.
2. Magnetic rods placed in back of brace have the purpose of support.
3. Hook and loop fasteners for easy application.

INDUSTRIAL BELT

07U IN - INDUSTRIAL BELT

Back part 24 cm / 9.5" high.

Aids in the prevention of back pain incurred during strenuous work.

Worn over clothing.

The detachable suspenders enable the person to wear the belt open.

Close the belt each time when lifting heavy objects.

In case of expected exertion, simply fasten the belt with the velcro fasteners.

Colors:

w

b

c



07U IN



07U - BACK BELT

Back part 24 cm / 9.5" high.

1. Useful for treatment of pain in the lower back.
2. Made of bi-elastic, mat.
3. Equipped with 4 support rods situated along the spine.
4. Belt alleviates tension on the spine by controlling abdominal sag.
5. 2 support rods along the spine.

07U




077A - BACK BELT

ES ITEM 07U

Back part 24 cm / 9.5" high.

Made of airy fabric

077A



078 - BACK BELT

ES ITEM 07U

Back part 20 cm / 8" high.

078





URIEL

Sport

HIS

PRESSURE
GARMENT

16



URIEL

BELT



07AU

07AU - LUMBO SACRAL (L-S)B.

Back part 24 cm 9.5" high.

1. Useful for treatment of pain in the lower back.
2. Made of bi-elastic, knitted cotton light pore material.
3. Equipped with 4 support rods situated along the spine.



T08

T08 - CRISS CROSS BELT

Back part 25 cm / 10" high.

Back belt with two crossed elastic straps on the back side for increased tension.

Useful for the treatment of mechanical low back pain.

Made of:

NEOPRENE LAMINATE



066

066 - LUMBO SACRAL BELT

1. Useful for treatment of pain in the lower back.
2. The 4 metal rods along the spine ensure maximum support.

NEW



0622

0622 - LUMBOSIL SACRAL BELT

1. Useful for treatment of pain in the lower back.
2. With silicone pelota.

Made of:

bi-elastic knitted material

NEW

064 - ATHLETIC BELT.

Back part 24 cm / 9.5" high.
Front part 16 cm / 6.5"

1. Useful for treatment of pain in the lower back.

with pelota



064



07A - LUMBO SACRAL (L-S)B.

Back part 24 cm / 9.5" high.

1. Useful for treatment of pain in the lower back.
2. Made of bi-elastic, knitted cotton light pores material.
3. Equipped with 4 support rods situated along the spine.

07A



07722 - LUMBO SACRAL (L-S)Belt

with SILICONE pelota

Back part 24 cm / 9.5" high.

1. Useful for treatment of pain in the lower back.
2. Made of bi-elastic, knitted Light porous material.
3. Equipped with 4 support rods situated along the spine.

0722




with silicone cushion

07AF - LORDOSIS BELT - FORTE

Back part 25 cm / 9.5" high.

1. Useful for treatment of strong pain in the lower back.
2. Made of bi-elastic knitted cotton, light porous material.
3. Straps adjust in the back and in the front.
4. Equipped with 8 support rods.
5. Provides extra strong support .

07AF



**PRESSURE
GARMENT**



18 

URIEL
SUPPORTO



FootCare

URIEL
Sport

HIS

**PRESSURE
GARMENT**



03

03 - SACRO-ILIAC-BELT

Supports the hips and stabilizes the sacroiliac joints. Recommended for pregnant women with problems caused by pregnancy.

Made of leather.



50

50 - SUSPENSORY

1. Supports swollen or sagging testicles.
2. Supports the testicle while the penis extends through an opening in the pouch.
3. Conditions requiring this support are injuries to the testicles, post surgical conditions, varicosity of the spermatic veins or an inherent weakness.
4. Light material designed for warmer climates.



59D

59D - HERNIA BELT - SUITABLE FOR BOTH SIDES

1. Suitable for inguinal hernia and post operative recuperation.
2. The belt includes 2 anatomically shaped foam pads which are removable thus allowing the belt to be used for the left + right side.
3. The belt has no metal parts (allows showering or swimming).
4. For post operative hernia recovery, please note that both foam pads must be removed.

70



70 - PAVLIK BANDAGE

For treatment of dysplasie / luxation
Available in 4 sizes

72



72 - ABDUCTION SUPPORT (mesh)

For treatment of dysplasie / luxation
Available in 6 sizes.

499 - HIP PROTECTOR

The Hip Protector can considerably reduce the risk of hip fracture.

The 2 shells – inner/outer shell - placed in short legged underpants, are designed to absorb the impact and spread the effect over a large area around the thighbone.

NEW



49

49 - THIGH SLEEVE T49 - THIGH SLEEVE

T49

1. Firm support for strained thigh muscles.
2. Suitable for all sporting activities.
3. Firm compression.
4. Protection against potential injuries, protective device.



MADE OF :
NEOPRENE LAMINATE

49PS - HIP BANDAGE

Post-operative bandage.
Keeps stitches and pads in place.

MADE OF :
LIGHT VENTILATED COTTONE MATERIAL

49 PS



49PSF - HIP BANDAGE FORTE

Post-operative bandage.

With spiral metal stays giving
extra support and stability

NEW

49PSF



**PRESSURE
GARMENT**



20



KNEE-PATELLA SUPPORT



T402

T402 - PATELLA FIX SUPPORT

1. Suitable for Patellafemoral condition.
2. Retropatellar arthrosis.

Suitable for both sides.

Made of: NEOPRENE



40

40 - PATELLA STRAP BANDAGE

4022 - PATELLA STRAP BANDAGE

with silicone cushion

T40 - PATELLA STRAP BANDAGE

Made of: NEOPRENE



with silicone cushion

4022

ONE SIZE FITS ALL

1. Suitable for Chondromalacia - patella condition.
2. Pressure directed to the patella tendon, support for a straight movement of the patella.
3. Relief for painful and fatigued knees.
4. Useful as support brace for minor knee complaints.



Made of: NEOPRENE

T40



T401

T401 - PATELLA STRAP KNEE SUPPORT

Combination of item T40 and T45X.

1. Suitable for Chondromalacia patella condition.
2. Pressure directed to the patella tendon, support for a straight movement of the patella.
3. Relief for painful and fatigued knees.
4. Useful as support brace for minor knee complaints.


Made of: NEOPRENE

FootCare

URIEL Sport

HIS

PRESSURE GARMENT



**PRESSURE
GARMENT**

41x - BICYCLE KNEE BRACE

41 - CRISS CROSS KNEE BRACE

1. Relief for painful, fatigued knees.
2. Support brace for minor knee complaints.
3. The "cross cut" feature helps to prevent irritation caused by flexing the knee. This brace exerts minimum pressure to the back of the knee.
4. 2 Spiral metal side stays for an extra support and stability.

(Length: 41 cm)



41X



45 - KNEE SLEEVE
Universal shape.

(Length: 24cm)

T45 - KNEE SLEEVE
(Length: 28 cm)

T45X - KNEE SLEEVE
As item T45 with open Patella.

MADE OF :
NEOPRENE LAMINATE



45E - LONG KNEE SLEEVE

As item 45
(Length: 28 cm)



48 - THERMO KNEE - SLEEVE

(Length: 24cm)



T45



T45X



T39



39

T39 - SHIN SLEEVE

Relieves pain along the medial or lateral side of the lower leg, caused by aerobic dancing, running, hiking and continued foot impact on hard surfaces.

MADE OF :
NEOPRENE LAMINATE

22


URIEL®
SILICONE KNEE SLEEVE


422

422 - GENUSIL KNEE SLEEVE

Low pressure silicone knee pad prevents / decreases swelling.

1. This support is equipped with a silicone pad that relieves pressure on the patella and helps to decrease the swelling around the patella.
2. Light lateral spirals help to stabilize the knee joint.

SILICONE PAD


422+

422+ - GENUSIL KNEE SLEEVE

As item 422 but with circumference strap, gives more side support.

SILICONE PAD


422D

422D - GENUSIL HINGED KNEE SLEEVE

1. This support is equipped with a silicone pad that relieves pressure on the patella and helps to decrease swelling around the patella.
2. Lateral hinges help to stabilize the knee joint.

SILICONE PAD


423

423 - LONG GENUSIL KNEE BRACE (length:40 cm / 16")

423 R for right
423 L for left

Circumferential elastic strap

Silicone pad

Open patella

2 spirals rods

This support is equipped with a silicone pad that relieves pressure on the patella and helps to decrease the swelling around the patella. Light duty lateral knee stabilization hinge.


FootCare
URIEL Sport
HS
PRESSURE GARMENT
MEASUREMENTS

Measurements on page 43



T43SKI - SKI BANDAGE

1. Most suited to prevent sport injuries and pain in the knee area.
2. Used in prevention of injuries (Neoprene impact-absorbing shell).
3. Provides increased heat retention.
4. Two spiral metal stays are inserted at each side of the knee giving extra support and stability.

new

T43 SKI



**T433D - HINGED KNEE SUPPORT
- DUAL CENTRIC HINGED**

(Length: 34 cm)

1. With fixed D hinges (non removable) to stabilize the knee joint. Useful in cases of meniscus and knee injuries.
2. Supplied with an adjustable horse-shoe shaped cushion.

T433D



All the supports supplied with adjustable horse-shoe shaped cushion which is fitted on, to help hold the patella in right position.



**PRESSURE
GARMENT**



24



HINGED KNEE BRACE



FootCare

URIEL
Sport

HIS

PRESSURE
GARMENT



43D

43D - HINGED KNEE BRACE (Length: 28 cm/11.2")

Suitable for persons who are unable to bend down / bend the leg.

With hook and loop fasteners in front of brace for easy application.

1. Helps restrict medial and lateral knee movements.
2. Supports and protects injured ligaments.
3. Useful for treatment of cases of arthritis and as a post operative brace.
4. Anatomical design, hinged knee joint.



T43VD

T43VD - HINGED KNEE BRACE (Length: 34 cm/ 13.5")

Suitable for persons who are unable to bend down / bend the leg.

With hook and loop fasteners in front of brace for easy application.

1. Helps restrict medial and lateral knee movements.
2. Supports and protects injured ligaments.
3. Useful for treatment of cases of arthritis and as a post operative brace.
4. Anatomical design, hinged knee joint.



T44G

T44G - KNEE SPLINT (Length: 40 cm/ 16")

HINGED KNEE BRACE
WITH FRONT
CLOSURE (40 cm)
ADJUSTABLE HINGES

1. Used for the support of the knee after operation, it serves as an alternative to a plaster-cast.
2. Equipped with adjustable anatomically shaped hinges at the kneejoint.

44 - KNEE IMMOBILIZER

(Length: 57 cm)

Used for the support and immobilization of the knee after operation. It serves as an alternative to a plaster-cast. Made of light material, enables ambulation. Lining made of terry cotton. Necessary treatments (removal of stitches etc.) can be easily continued. The hook & loop fasteners make it possible to fit the bandage according to the various, desired adjustments. The posterior rods are preshaped & provide rigid immobilization. The bandage can be washed (after removal of rods).




available in jeans-color

44G - KNEE SPLINT DUALCENTRIC HINGED

with adjustable extension / flexion hinges
(Length: 53cm **S-M** 57cm **L-XL** 57cm **XXL**)




1. Used for the support of the knee after operation, it serves as an alternative to a plaster-cast.
2. Equipped with adjustable anatomically shaped hinges at the knee joint.
3. Enables extension, flexion and fixation.
4. Made of light material.
5. Necessary treatments (removal of stitches etc.) can be continued easily.
6. The hook and loop fasteners make it possible to fit the splint according to the various, desired adjustments.



DUALCENTRIC
G HINGE

RECOMMENDED TO USE DUALCENTRIC D/G HINGE





SYSTEM OF
DUALCENTRIC
HINGE D/G



HING G

DUALCENTRIC
HING +
ADJUSTABLE
EXTENSION AND
FLEXION
HINGE

SYSTEM OF
SINGLE AXIS
HINGE H



FootCare

URIEL
Sport
HIS

PRESSURE
GARMENT


**PRESSURE
GARMENT**
30 - ANKLE SUPPORT (adjustable)
30

1. Helps ease painful, tired feet and ankles.
2. Useful as support for strained and sprained ligaments.
3. Hook and loop fasteners for easy application.
4. Adjustable sizes and tension due to fasteners.
5. Applicable to left or right foot.


35 - ANKLE SUPPORT
35 C

1. Helps ease painful, tired feet and ankles.
2. Useful as support for strained and sprained ligaments.
3. Applicable to left or right foot.
4. Extended length of the support.
5. Useful in minor cases of ankle injuries.


35

34 - ANKLE SOCK
34

1. Helps ease painful, tired feet and swollen ankles.
2. Useful as support for strained and sprained ligaments.
3. Applicable to left or right foot.
4. Covers the heel totally.
5. Extended length of the support.


34+


As item 34 with additional
fixation strap for extra support.



27



ANKLE SUPPORT



silicone cushion

322 R

322 L

ANKLE SILICONE BANDAGE

322 R - for right side

322 L - for left side

Ankle support, for treatment of lateral swelling of the malleo bone and postoperative swelling. With silicone cushion and fastener closure for correct adjustment at any stage of the swelling.

with silicone cushion.



silicone cushion

HTS 322 R

HTS 322 L

MALLEOFLEX SUPPORT

HTS 322R - for right

HTS 322L - for left

Ankle support for the treatment of lateral swelling of the Maleo bone.

With silicone cushion and fastener closure for correct adjustment at any stage of the swelling.

Made of thin Lycra fabric, requires minimum space and therefore fits in any shoe.




HTS 355

HTS 355 - ACHILLOFLEX SUPPORT

For relief of pressure on the achilles tendon. Comfort for a shortened achilles tendon. Rubber padding to avoid pressu on the Achilles sinew.

Extra heel cushion for other foot to prevent pelvic disalignment.

Made of thin Lycra fabric, requires minimum space and therefore suitable for any shoe.



333

333 - ANKLE AIR BRACE

Universal size. Fits right or left

The Air Brace provides functional ankle support, suitable for ankle injuries such as sprains and unstable ankles. The anatomical shaped plastic shell and the foamed air-bubble pad, worn with a laced shoe, effectively limits inversion/eversion, however normal flexion during walking and sporting activities is virtually unlimited. Furthermore it prevents mechanical ankle pains and ankle injuries.



**PRESSURE
GARMENT**

MEASUREMENTS

Measurements on page 43

33 - LIGHT ANKLE SPLINT

1. This splint stabilizes the ankle joint with two lateral plastic stays that have an open space over the Malleoli to prevent undue pressure.
2. Useful for the treatment of inflammation of the ankle joint and in cases of Tendinitis, Arthritis, strained tendons and sensitive ankle joints.



33

PERFECT FITTING

Made of:
STRONG KNITTED COTTON
ELASTIC MATERIAL

**COLOR
B**
T33


MADE OF NEOPRENE


**W
COLOR**

**C
COLOR**
INTSIDE

OUTSIDE

32


T32


MADE OF NEOPRENE

32 - ANKLE SPLINT

Conditions:

1. **Ideal for Basketball players** and other sports. Helps prevent, injuries to the foot/ankle joint.
2. After fractures. After removal of the plaster-cast, the splint provides side support.
3. Useful for the treatment of inflammation of the ankle joint and in cases of Tendinitis, Arthritis, strained tendons and sensitive ankle joints.
4. To stabilize the ankle after operation. The removable splint makes it possible to continue the necessary treatment.

This splint stabilizes the ankle-joint with 2 supporting rigid splints on the medial side and 2 flexible spiral rods on the lateral side, without limiting the movement of the foot. Permits free activities.



**PRESSURE
GARMENT**



29



WALKERS



URIEL
Sport



PRESSURE
GARMENT

MEASUREMENTS

1300

1300 - POST-OP-SHOE

Deal for use following surgery, as well as for post trauma, dislocations and fractures in the forefoot area. The hook and loop closure provides even tension and compression. The rigid outsole provides stability and minimizes flexure, soft breathable upper. Fits right or left.

1303

1303 - FOREFOOT RELIEVER

Designed for forefoot relief and healing after surgery of the toes. This shoe has a wedge outsole and a soft breathable upper. The hook and loop closure provides even tension and compression. Fits right or left.

1305

1305 - CAST SANDAL

Fits on walking cast. Provides walking comfort.

1395

1395 - WALKING CAST

Provides an alternative to traditional casting methods, for the treatment of fractures of the lower leg, foot and ankle and ankle sprains. May be adjusted to accommodate for changes due to swelling, which permits consistency if immobilization. Designed for maximum patient comfort, permits examination, skin care and proper hygiene. With foam leg-wraps .

1396

1396 - WALKING CAST

With adjustable Hinge

As item 1395.
May be adjusted for varying degrees of plantar or dorsi-flexion, for controlled rehabilitation.



Adjustable Hinge



30



URIEL



URIEL

URIEL
Sport
HISPRESSURE
GARMENT

MEASUREMENTS

**364 S-M-L**

1 Pair

364 XL

1 Piece

TOE SPREADER

Silicone gel pad, which contains mineral oil. Avoids the overlapping of toes, prevents, reduces and soothes pain caused by calluses.

**312 9X5 cm, 4 mm thickness
ADHESIVE SILICONE SHEET**

1 Piece

Adhesive silicone sheet for toes / feet, which contains softening mineral oil. It prevents pressure and rubbing on calluses and sensitive skin. Guarantees quick pain relief.

**366****CORN PAD**

2 Piece

Adhesive silicone pad for toes/feet. It prevents pressure and rubbing on calluses and sensitive skin. Guarantees quick pain relief.

**3661****CORN REMOVER**

6 Piece

Active plaster for the removal of corns. Adhesive backing keeps the disk in position, the inserted protective ring avoids pressure and friction and provides perfect comfort. The disk contains salicylic acid for the treatment of callosity, the ingredients soften the corn or callus and permit their removal.

**365 S-M-L**

1 Pair

365 XL

1 Piece

CORN PAD PROTECTOR

Elastic ring with silicone gel corn-pad, which contains mineral oil. It prevents, reduces and soothes pain caused by rubbing. Recommended for blisters, calluses, corns, hammer-toes, etc.

**3655 S-M-L
TOE TUBE**

1 Piece

Tube coated with silicone gel provides immediate cushioning, for toes and toe nails prevents friction and pressure .not to be used on open wounds, washable and reusable ,cut to desired length. contents: 1 tube, 15 cm length

**369 S-M-L**

1 Piece

369 XL

1 Piece

TOE/FINGER TUBE CUSHION

Cap of elastic fabric lined with silicone gel, which contains mineral oil. Reduces, removes and prevents toe or nail lesions.

367 One size
BUNION CUSHION

1 Piece

Silicone gel pad, which contains mineral oil. It prevents, reduces and soothes pain caused by bursitis and bony exostitis of the first metatarsal-phalanx joint.



3677
BUNION PROTECTION RING

2 Pieces

Adhesive silicone pad for the bunion. Prevents pressure and rubbing on calluses and sensitive skin. Guarantees quick pain relief.



3641 XL
TOE SPREADER

1 Piece

Made of rigid silicone material. Straightens the big toe. It avoids the overlapping of toes. Prevents, reduces and soothes pain caused by calluses.



+

388R For right foot. 388L For left foot.

One size

HALUX VALGUS NIGHTTIME SPLINT

To prevent surgery of the halux valgus. Prevents abduction. It is advisable to wear art.366 during the day.



38 - HALUX VALGUS BANDAGE

Comfortable after surgery of the Halux Valgus. Prevents abduction.

Fits left and right. One size fits all.



311 One size
ACTIVE GEL PLASTER

5 Pieces

For the treatment of blisters. Stops friction and relieves pain. Highly conformable cushioning, breathable and waterproof plaster. The hydrocolloid dressing promotes rapid healing.





32



URIEL



URIEL

URIEL
Sport
HISPRESSURE
GARMENT

MEASUREMENTS

**3623 One size
TOE REGULATOR**
1 Piece

Foam cushion for the ball of the foot with elastic loop for the second toe. Protects hammer-toes from rubbing against the top of the shoe.

**368 One size
HAMMER-TOE CUSHION**
1 Piece

Special foam cushion. Prevents pressure on sensitive areas of the toe.

**363 One size
TOE CUSHION**
1 Pair. For right and left foot

Silicone gel pad, which contains mineral oil. It avoids calluses and hammer-toes and stretches contracted toes. With fastening ring for third toe.

**361 One size
FOREFOOT CUSHION BANDAGE**
1 Pair

Elastic bandage with silicone gel pad, which contains mineral oil. Guarantees a shock absorbing action. Protects the forefoot area (metatarsals) from pressure and friction.

**362 One size
METATARSAL CUSHION**
1 Pair

Silicone gel pad, which contains mineral oil. It protects ligaments, tendons and joints. Avoids the formation of calluses and prevents inflammation in the forefoot.

**3622 One size
CALLUS PAD**
4 Piece

Adhesive silicone pad for the feet, which contains softening mineral oil. It prevents pressure and rubbing on calluses and sensitive skin. Guarantees quick pain relief.



305W / 305M For shoe sizes 35-46
IMPROVED ARCH SUPPORT
1 Pair

The ideal support for sports and daily use. For relief of pain in foot and lower back. For equal distribution of body weight.

Improved insole with 2 removable P.V.C.-cushions of different thickness, supports the inner arch due to individual request.

NEW



377 For shoe sizes 35-46
ORTHOPEDIC INSOLE
1 Pair

Orthopedically shaped insole made of latex, coated on one side with 100% cotton towelling. Relieves pain in foot and back area, especially suitable for sport-shoes.

1 Pair



303 For shoe sizes 36-45
SILICONE DIABETIC ARCH SUPPORT WITH OUT M.T. PAD

Soft cushioning. Protection of sensitive soles of diabetics. Pressure relief.


304 For shoe sizes 36-45

1 Pair

SILICONE DIABETIC ARCH SUPPORT

WITH M.T. PAD + REPLACABLE ARCH PADS
Soft cushioning. Protection of sensitive soles of diabetics. Pressure relief.

REPLACABLE ARCH PADS



METATARSAL SUPPORT

360 One size fits all

1 Pair

SILICONE FOREFOOT RELEASER

Narrow and light spread-foot support with integrated metatarsal pad. Requires minimum space, suitable for any shoe. The metatarsal pad provides an anatomic support for the forefoot arch and prevents pressure in the forefoot.



REMOVEBEL M.T. PAD



306 For shoe sizes 36-45

1 Pair

SILICONE HEEL CUSHIONS

307 For shoe sizes 36-45

1 Pair

SILICONE HEEL CUSHIONS

Soft blue dot in the center to relieve pain in the heel bone area, reduced pressure on the spur during walking.



309 For shoe sizes 36-45

1 Pair

SILICONE HEEL SUPPORT CUSHIONS

as item 307 with additional support of the inner arch.

NEW





34



URIEL

FootCare

URIEL



FootCare

URIEL
Sport
HISPRESSURE
GARMENT

MEASUREMENTS

**No FC 310 PROTECTIVE FOOT CREAM**

A natural preventive foot care ointment for maintaining hygienic conditions around the feet. Contains ingredients such as tea-tree-oil, salicylic acid and Dead-Sea-Minerals. Ideal for delicate areas such as toes, for everyday use, especially in a warm or humid climate. Leaves the feet fresh and smooth. Recommended for athletes or before physical activity.

No FC 311 DEO FRESH FOOT CREAM

A deodorizing foot cream suitable for daily use. Antiperspirant action keeps skin dry and dispels undesired odours. The formula is enriched with an astringent extract, Dead-Sea- Minerals and aromatic oil derivatives for a sensation of freshness during the whole day.

No FC 312 SOFTENING FOOT SCRUB

Special cream for scrubbing and softening dry and rough foot-skin.

The exfoliating scrub contains natural tea-tree-grains and is enriched with aloe vera, plant extracts and Dead-Sea-Minerals for fresh and smooth feet.

No FC 313 SMOOTHING FOOT BALSAM

The herbal protective foot cream for rough, dry and cracked skin. Enriched with natural plant extracts selected essential oils and Dead-Sea-Minerals to refresh the skin, enhance it's softness and elasticity and maintain hygienic conditions.

No FC 314 COOLING LIQUID TALC

A refreshing foot care lotion which rapidly dries on the feet leaving a velvety soft foot powder. Provides a deodorizing effect with a pleasant cooling sensation. The lotion maintains the hygienic conditions of the feet leaving the skin dry and perfectly smooth.

No FC 315 ICY MINT COOLING GEL

A special cooling and calming gel that serves as an instant relief for tired legs. Contains a handful of herbal astringent such as extract of Vine Leaves; with Hazel, Cypress and Ruscus roots. All are known to have a toning effect on the capillaries and soothe congested extremities. The formula gives a refreshing sensation even after a single use. The gel may be used as often as required and is suitable for all kind of skin types.

No FC 330 SHOE DEODORANT

The shoe refresher allows you to spray inside all types of footwear quickly and easily, especially recommended for sport-shoes, army-boots and slippers. Eliminates and prevents unpleasant odours according to the new effective, anti-bacterial formula. Leaves footwear smelling clean and fresh. Works best if used immediately after removing shoes.



No FC 331 FOOT SPRAY

Especially efficient properties for active foot care. Contains antibacterial and antifungal ingredients which guarantee long-term protection. Clinically proven, antiperspirant action keeps the feet fresh and dry all day long.



No FC 341 NAIL FUNGUS DROPS

Liquid for the treatment of fungal nail conditions. After removing the infected parts apply the liquid on the nail.



No FC 342 INGROWN TOE NAIL DROPS

Softening imbedded toe nails and calloused tissue, prevents the ingrowing of the toenail and enables the straightening of the nails.

No FC 321 BATH SALTS

Foot-Bath Salts with Dead-Sea-Minerals. The solution accelerates the metabolism and increases the secretion of the skin. Cools and refreshes hot and aching feet. For relaxed feet and legs and for smooth and supple skin.



No FC 350 ANTISEPTIC FOOT-SOAP

Cleans gently and prevents foot-perspiration by the deodorizing and antiseptic effect of the Australian Teatree Oil and maintaining this condition for several hours. with Dead-Sea-Minerals. Refreshing and relaxing tired feet.



HTS
PRESSURE
GARMENT

MEASUREMENTS

36



S - SOCKS[®]

Silver socks



390 DIABETIC SILVER SOCKS

Pressure-free and super-soft special socks for diabetics. Provide maximum wearing comfort by seamless knitted fabric and best possible material composition. The incorporated silver has an antiseptic and deodorizing effect and compensates for changes in temperature. Perfect fitting, made without any pressing elastic.



391 ANTISEPTIC SILVER SOCKS

The socks create a permanent feeling of well-being for stressed feet. They are pleasant to wear and the incorporated silver has antiseptic and deodorizing properties. Compensates for changes in temperature. Specially suitable for people suffering from skin irritation.



392 SPORT SILVER SOCKS

The socks create a permanent feeling of well-being for stressed feet. They are pleasant to wear and the incorporated silver has antiseptic and deodorizing properties. Compensates for changes in temperature. Specially suitable for sports.

All socks available in colours:



**PRESSURE
GARMENT**

MEASUREMENTS

PEDICURE PRODUCTS



PD 340
Nail clipper



PD 345
Corn EX



PD 341
Corn cutter



PD 346
10 blades
for corn cutter



PD 342
Flat callus rasp
(coarse & fine)



PD 343
Callus rasp with
natural pumice stone

PD 344
Pumice stone



347
Nail nippers
for ingrown
toe-nails



348
Nail nippers




349
Nail nippers

INSOLES



370
Latex
Cushion

371
Terry
Comfort



372
Odour
Destroyer

374
Shock
Absorber

375
Thermo
Insole



**PRESSURE
GARMENT**



SPECIAL DEPARTMENT FOR THE PRODUCTION OF PRESSURE GARMENT

In this department we produce - from special elastic fabric and due to individual sizes:

1. Elastic gloves after hand dissections in order to accelerate the recuperation of the scars.
2. A wide collection of elastic pressure garment for use after burn injuries.

1. Pressure Garment - Burn Injuries



2. Glove - After Hand Dissection



Children - Pressure Garment - Collection



For the receipt of a product
catalogue kindly contact us.



PRESSURE
GARMENT

SPECIAL CORSETS, AFTER PLASTIC SURGERY LIKE
LIPOSUCTION



Measures before operation Maße vor Operation
Mesures avant l'opération Maß inden operation


in cm.

	XS	S	M	L	XL	XXL
	66 - 75	76 - 85	86 - 95	96 - 105	106 - 115	116 - 125
	60 - 69	70 - 79	80 - 89	90 - 99	100 - 109	110 - 119
	90 - 99	100 - 109	110 - 119	120 - 129	130 - 139	140 - 149

in inch.

	XS	S	M	L	XL	XXL
	26.5 - 30	30.5 - 34	34.5 - 38	38.5 - 42	48.5 - 46	46.5 - 50
	24 - 27.5	28 - 31.5	32 - 35.5	36 - 39.5	40 - 43.5	44 - 47.5
	36 - 39.5	40 - 43.5	44 - 47.5	48 - 51.5	52 - 55.5	47.5 - 59.5

LIPOSUCTION



**PRESSURE
GARMENT**

MEASUREMENTS

HI-TEX-SUPPORT IS MADE OF AN INNOVATIVE,
ESPECIALLY THIN LYCRA FABRIC.
OPTIMAL FITTING THE INDIVIDUAL ANATOMICAL SHAPE OF THE BODY LIKE A
SECOND SKIN ACCORDING TO THE MEMORY EFFECT

MATERIAL COMPOSITION:

42% LYCRA
16% POLYURITHAN
42% POLYMICRO
FLEECE



BREATHABLE
LYCRA FABRIC
POLYURITHAN
POLY FLEECE
BODY TEMPERATURE

•RETAINING BODY HEAT•BREATHABLE•PREVENTING FOLDS•



FootCare

URTEL
Sport
HTS

PRESSURE
GARMENT

Neoprene
bandage
support

URIEL^R
Sport



S-XL



41

ONE SIZE FITS ALL
UNE TAILLE POUR TOUS
EINE GRÖSSE - PASST ALLEN
TALLE PARA TODAS LAS MEDIDAS

SPECIAL RANGE FOR PHARMACIES



FootCare

URIEL^R
Sport

HIS

PRESSURE
GARMENT

MEASUREMENTS

URIEL^R
Sport





42

Item 60 - 66



S		32-37
M		38-43
L		44-48



Size	Item 80,82
S	CHILD
M	WOMAN
L	MAN

Size	Item 85,85+
S	CHILD
M	CHILD
L	WOMAN
XL	MAN



Item 90,91, 922,95,T98	In.	Size	Cm.
	8 - 9	XS	21-23
	9 - 10	S	23-25
	10-11	M	25-28
	11-12	L	28-30
	12-13	XL	30-32
	13-14	XXL	32-34

Item 78,86	In.	Size	Cm.
	8.0 - 9.2	XS	21-23
	9.2 - 10.8	S	23-25
	10.8 - 12.0	M	25-28
	12.0 - 13.8	L	28-30
	13.2 - 14.4	XL	30-32
	14.4 - 16.0	XXL	32-34



Item 20 - 28	In.	Size	Cm.
	4.5-5.5	XS	12-14
	5.5-6.5	S	14-16
	6.5-7.5	M	16-18
	7.5-8.5	L	18-20
	8.5-9.5	XL	20-22
	9.5-10.5	XXL	22-24



Item 74,75,76	In.	Size	Cm.
	10-24	XS	50-60
	24-28	S	60-70
	28-32	M	70-80
	32-36	L	80-90
	36-40	XL	90-100
	40-44	XXL	100-110

Item 776,79,87,	In.	Size	Cm.
	26-30	XS	65-75
	30-34	S	75-85
	34-38	M	85-95
	38-42	L	95-105
	42-46	XL	105-115
	46-50	XXL	115-125



Item 06-12	In.	Size	Cm.
	26-30	XS	65-75
	30-34	S	75-85
	34-38	M	85-95
	38-42	L	95-105
	42-46	XL	105-115
	46-50	XXL	115-125

Maternity belt		Item 18	
Clothes Size	U.S.A	Size	Europe
	4-6	S	36-38
	8-12	M	40-42
	10-12	L	44-46
	14-16	XL	48-50
	18-20	XXL	52-54

Maternity belt		Item 19	
Clothes Size	U.S.A	Size	Europe
	4-10	S -M	36-40
	10-16	L-XL	42-44
	18-20	XXL	46-48

DRIEL
Sport

HIS


PRESSURE
GARMENT

Item 70	Size	In.	Cm.
	S	13.5-16.0	34.0-40.0
	M	16.5-19.5	16.5-19.5
	L	20.0-23.5	20.0-23.5

Item 72	Size	In.	Cm.
	XS	4	10
	S	6	15
	M	8	20
	L	10	25
	XL	12	30
	XXL	14	35

Item 03, 499,49PSF,49PS 50, 55, 59	Size	In.	Cm.
	XS	26-30	65-75
	S	30-34	75-85
	M	34-38	85-95
	L	38-42	95-105
	XL	42-46	105-115
	XXL	46-50	115-125


Item 49	Size	In.	Cm.
	XS	12-15	30-37
	S	15-18	37-45
	M	18-21	45-52
	L	21-24	52-60
	XL	24-27	60-68
	XXL	27-30	68-75


Item 41 - 48	Size	In.	Cm.
	XS	8.5-10.5	22-27
	S	10.5-12.5	27-32
	M	12.5-14.5	32-37
	L	14.5-16.5	37-42
	XL	16.5-18.5	42-47
	XXL	18.5-20.5	47-52

Item 44	Size	In.	Cm.
	S-M	15.5-17.5	39-44
	L-XL	18.0-22.5	45-57
	XXL	23.0-24.5	56-62

Item 30,32,33,35	Size	In.	Cm.
	XS	6-7.25	17-19
	S	7.25-8	19-21
	M	8-9.25	21-23
	L	9.25-10	23-25
	XL	10-11.25	25-27
	XXL	11.25-12	27-28

Item 322,34,34+, 35	Size	In.	Cm.
	S	7.25-8	19-21
	M	8-9.25	21-23
	L	9.25-10	23-25
	XL	10-11.25	25-27





ITEM 1012
69x97cm



ITEM 1005
34x97cm



ITEM 1004
45x60cm

POSTERS



ITEM 1000 U
51X30X200



PHARMACY DISPLAY



ITEM 1003 FC
23X18x62cm

FOOT-CARE DISPLAY



Foot-care turnabel Display
ITEM 1002 FC



Foot-care Display Stand
ITEM 1000 FC



PRESSURE
GARMENT



FootCare

DRIEL
Sport
HIS


PRESSURE
GARMENT

MEASUREMENTS



SCENES OF MACHINERY FROM THE FACTORY


SCENES OF MACHINERY FROM THE FACTORY



Circular tube knitting machine



CAD-computer design system





48



URIEL
FIRST AID



FootCare

URIEL
Sport

HIS

PRESSURE
GARMENT

MEASUREMENTS

ITEM 801



ICE & GO

ELASTIC COLD BANDAGE

NEW!

- EASY TO USE
- COOLS FOR HOURS
- MAXIMUM COMFORT
- IMMEDIATE RELIEF

URIEL
FIRST AID